

OpenSignals (r)evolution File Formats Description

OSF 080517

FEATURES


- > Data visualization & recording
- > Feature extraction plugins with reporting
- > Multiplatform support
- > User-friendly GUI

APPLICATIONS

- > Psychophysiology
- > Biomedical projects
- > Computer science
- > Electrical engineering
- > Human-Computer Interaction
- > Robotics & Cybernetics
- > Physiology studies
- > Biomechanics
- > Biofeedback

GENERAL DESCRIPTION

OpenSignals records data in ASCII tab delimited files (Fig. 1) and/or in the modern and highly efficient HDF5 format (Fig. 2). If the same format is followed, any third-party software can also export or record data in such a way that OpenSignals may be used as visualizer.



```
# OpenSignals Text File Format
# ("20:16:04:12:01:93": {"sensor": ["RAW", "RAW", "RAW", "RAW", "RAW", "RAW"], "device
name": "20:16:04:12:01:93", "column": ["nSeq", "I1", "I2", "O1", "O2", "A1", "A2", "A3",
"A4", "A5", "A6"], "sync interval": 2, "time": "20:56:46.178", "comments": "", "device
connection": "20:16:04:12:01:93", "channels": [1, 2, 3, 4, 5, 6], "date": "2017-4-18",
"model": 0, "digital I/O": [0, 0, 1, 1], "firmware version": "5.2", "device":
"bitalino_rev", "position": 0, "sampling rate": 1000, "label": ["A1", "A2", "A3", "A4",
"A5", "A6"], "resolution": [4, 1, 1, 1, 1, 10, 10, 10, 10, 6, 6], "special": [{"}, {}, {}],
{}}, {}})
# EndOfHeader
1 1 1 0 0 511 495 0 534 38 4
2 1 1 0 0 506 495 0 529 38 5
3 1 1 0 0 516 496 0 527 38 5
4 1 1 0 0 514 498 0 535 38 6
5 1 1 0 0 515 499 0 530 38 6
6 1 1 0 0 511 500 0 569 38 5
7 1 1 0 0 499 499 0 579 38 5
8 1 1 0 0 515 500 0 576 38 3
9 1 1 0 0 515 502 0 571 38 2
10 1 1 0 0 511 504 0 575 38 3
11 1 1 0 0 506 507 0 579 38 4
12 1 1 0 0 516 507 0 577 38 5
13 1 1 0 0 515 506 0 576 38 5
14 1 1 0 0 514 506 0 572 38 6
15 1 1 0 0 511 507 0 566 38 6
0 1 1 0 0 521 507 0 578 38 5
1 1 1 0 0 508 507 0 573 38 5
2 1 1 0 0 515 507 0 568 38 3
3 1 1 0 0 511 503 0 563 38 2
4 1 1 0 0 513 502 0 561 38 3
```

Fig. 1. Example ASCII file for BITalino (r)evolution; **Col. 1:** Sample sequence number (4-bit) generated on the device to facilitate the detection of missing data; **Col. 2 & 3:** State of the digital inputs I1 & I2; **Col. 3 & 4:** State of the digital outputs O1 & O2; **Col. 5-:** Raw data for each analog input acquired during the recording session (in this case A1-A6).




Fig. 2. Example of a recording stored in HDF5 format.



PLUX – Wireless Biosignals, S.A.
Av. 5 de Outubro, n. 70 – 8.
1050-059 Lisbon, Portugal
plux@plux.info
<http://www.plux.info/>

REV A

© 2017 PLUX 

This information is provided "as is," and we make no express or implied warranties whatsoever with respect to functionality, operability, use, fitness for a particular purpose, or infringement of rights. We expressly disclaim any liability whatsoever for any direct, indirect, consequential, incidental or special damages, including, without limitation, lost revenues, lost profits, losses resulting from business interruption or loss of data, regardless of the form of action or legal theory under which the liability may be asserted, even if advised of the possibility of such damages.

OpenSignals (r)evolution

File Formats Description

OSF 080517

METADATA FIELDS SUMMARY

Field	Description
<i>channels</i>	Set of analog inputs selected for acquisition in the recording session that generated this file (last columns of the file)
<i>column</i>	Meaning of each column (in ASCII files) or group with acquired data (in HDF5 files) for a given device ¹
<i>comments</i>	User-defined text with notes or other content of interest to the file
<i>date</i>	Day, month and year in which the file was recorded
<i>device</i>	Type of recording device used to collect the data (e.g. biosignalsplux, BITalino (r)evolution)
<i>device connection</i>	Logical address used to establish a connection with the device ²
<i>device name</i>	Friendly name manually assigned by the user to the device
<i>digital IO</i>	Set of digital channels available in this device (0 – Input; 1 – Output)
<i>firmware</i>	Version of the embedded software running on the device that acquired the data
<i>label</i>	Array with labels manually assigned by the user to each of the analogs channels
<i>mode</i>	Specifies the acquisition mode in which the device was used (0 – Regular acquisition; 1 – Multiple synchronized devices; 2 – Started by a triggering signal)
<i>resolution</i>	Resolution with which each of the columns (in ASCII files) or group with acquired data (in HDF5 files) are represented
<i>sampling rate</i>	Number of samples per second with which data has been acquired by the device on each channel
<i>sensor</i>	Array with the type of sensor connected to each analog port (e.g. as selected from the options available in the device configuration panel on OpenSignals)
<i>special</i>	List of non-standard channels acquired (e.g. SpO2 sensors with I2C interface connected on biosignalsplux hubs).
<i>sync interval</i>	Time interval (in seconds) at which a digital signal is sent by a “pacemaker” thread to a single device (used when the sync mode in on OpenSignals for synchronized data acquisition using multiple devices)
<i>position</i>	Order in which the block of columns corresponding to a given device appears in each line of the file (0 – First sequence of columns; 1 – Second sequence of columns; etc.)
<i>time</i>	Time at which the first sample was received by the software

ASCII TEXT FORMAT

OpenSignals ASCII text files have two parts, namely a header section and data section.

The header section has 3 header lines, each starting with the pound (#) character, given that it is automatically ignored by ASCII text file loading functions in common scientific computing programming languages (e.g. `loadtxt` in Python).

The header section begins after the line with the content “# OpenSignals Text File Format” and ends after the line “# EndOfHeader”. The second line has the metadata stored as a JSON object containing each device's MAC address as top-level keys. For each top-level key, there is a JSON object as value, containing the acquisition settings used for that

¹ i.e. **nSeq** is the sample sequence number generated on the device, **Ix** is the state of digital input x, **Ox** is the state of digital output x, and **Ax** is the raw data sampled by the device.

² e.g. **AA:AA:AA:AA:AA:AA** represents a connection over a Bluetooth socket, **0.0.0.0:0000** represents a connection over TCP/IP, and **COM1** or **/dev/tty.BITalino-AA-AA-DevB** represents a connection over serial port.


OpenSignals (r)evolution

File Formats Description

OSF 080517

device in the recording session that produced the file in the form of key / value fields (described earlier).

JSON has been chosen as the metadata representation format, given that in most scientific computing languages it can be directly unraveled as a dictionary, and it has parsers available for most mainstream programming languages.



```
# OpenSignals Text File Format
# {"20:16:04:12:01:93": {"sensor": ["RAW", "RAW", "RAW", "RAW", "RAW", "RAW"], "device
name": "20:16:04:12:01:93", "column": ["nSeq", "I1", "I2", "O1", "O2", "A1", "A2", "A3",
"A4", "A5", "A6"], "sync interval": 2, "time": "20:56:46.178", "comments": "", "device
connection": "/dev/tty.BITalino-01-93-DevB", "channels": [1, 2, 3, 4, 5, 6], "date":
"2017-4-18", "mode": 0, "digital I0": [0, 0, 1, 1], "firmware version": "5.2", "device":
"bitalino_rev", "position": 0, "sampling rate": 1000, "label": ["A1", "A2", "A3", "A4",
"A5", "A6"], "resolution": [4, 1, 1, 1, 1, 10, 10, 10, 10, 6, 6], "special": [{}, {}, {}
{}], {}, {}}}}
# EndOfHeader
```

Fig. 3. Example of an ASCII text file header.

On the data section, each column contains the signals labeled in the `column` field (see the caption on Fig. 1 for an example). To facilitate automated post-processing of the data, the resolution with which each column is recorded is stored in the `resolution` field (e.g. useful for a specific sensor's transfer function found in the corresponding sensor datasheet).

If the acquisition includes multiple devices, the metadata is stored following the same logic for each device, side by side, according to the `position` header field.

HDF5 FORMAT

OpenSignals HDF5 files comprise one group for each device used, represented by its MAC address. Each device's group includes metadata fields containing the acquisition settings used for that device in the recording session that produced the file (described earlier) and 5 sub-groups:

- > **digital**: includes a dataset for each digital channel (input channels followed by the output channels);
- > **events**: includes a dataset for digital events and dataset for sync events;
- > **plugin**: includes a group of datasets for each plugin used/processed;
- > **raw**: includes a dataset for the sample sequence number generated by the device and a dataset for each analog input selected for acquisition (each analog dataset includes the label and sensor type attributes).
- > **support**: includes a group of datasets with support information (mean, standard deviation, ...) for each zoom level available for each channel, either analog or digital (`t`: initial time for the sample group; `mean`: average, `sd`: standard deviation; `mx`: minimum value; `Mx`: maximum value; `mean_x2`: average of the sample group's 2nd power).

OpenSignals (r)evolution File Formats Description

OSF 080517




Fig. 4. Example of an HDF5 file structure.